Soft Band Patient Identification

A better patient ID solution

The *DT SoftBand* features the ability to print wristbands utilising your existing thermal laser printers.

Efficient. Safe. Clinically Effective.

Compared to other solutions, the direct thermal printer friendly **DT SoftBand** speeds the hospital admission process, whilst providing a smooth, comfortable identification solution for all users – from premature babies to elderly patients. With our patented phthalate and latex free material, our bands provide smear & smudge-free usage ensuring information is always easily read.

Customised to your Requirements.

The **DT SoftBand** is available in different sizes and styles for every patient you treat. From newborns to bariatric patients, we have the solution for the most complex organisation. Extenders can be supplied to give you a wider range of fittings.

Compliant with all Safety & Quality Specifications Worldwide.

The *DT SoftBand* patient identification wristbands can assist in patient identification to ensure every patient is safely and properly identified through the hospital experience. It allows instant identification of the patient direct from your Patient Information System.

Uses Existing Thermal Direct Printers.

DT SoftBand has been tested, certified and approved on a wide range of leading brand direct thermal printers, so hospitals are not locked into a single proprietary system.





For more information or for a sample, contact

Telephone: 1300 813 931

Email: info@designedtostick.com.au

www.designedtostick.com.au



SoftBand Direct Thermal Comfort Wristbands

Healthcare providers know the importance of keeping their patient comfortable. Our soft direct thermal wristbands offers a sleek shape and soft smooth material. We listened to customer requests for increased image space to fit critical patient identification. The band fastens securely with an aggressive self adhesive closure and tamper evident feature. The wristband is engineered to stay on during the entire stay of the patient. The image space is protected by a special coating to resist alcohol, hand sanitizer, and soapy water. Our product allows use of any common direct thermal printer so hospitals are not locked into a single proprietary system designed to prevent competition.









Ease of use. Simply peel the band off the carrier from the non-adhesive end and affix it onto the patient with the self adhesive closure. Patients will hardly notice they are wearing a wristband due to it's light weight comfortable material. Operations will appreciate the great pricing of this patient friendly band. Sheets of extenders individually separable by perforations are available to eliminate daisy chains of wristbands.



Available in Adult; Pediatric, and Baby sizes as shown or in custom shapes and sizes.









For more information or for a sample, contact

Telephone: 1300 813 931

info@designedtostick.com.au

www.designedtostick.com.au

